

Lucerna DR Series

Floor & Ceiling mounted DR



Lucerna DR 65
Ceiling



The new X-RAY generation
Digital Radiography Solution

Lucerna DR 65 Ceiling

VERY VERSATILE

Totally motorized, it allows both automatic and manual movements, providing high flexibility in use. The functions of autotracking and stitching make **Lucerna DR** a complete system for any type of examination. The possibility of using a double FPD makes it extremely versatile.



● Automated tilting Bucky



ACQUISITION SOFTWARE

- Management of the entire workflow.
- Complete control of digital images.
- Outstanding image quality.
- Automatic APR setting .
- Fully automatic Stitching.
- DICOM 3.0.

Lucerna DR 65
Floor Mounted



The new X-RAY generation
Digital Radiography Solution

Lucerna DR

Digital DR System

LUCERNA DR65F "Soft Handling" system, through the use of the "handle bar", with easy movements management. Equipped with a backlit inclinometer, it allows checking the inclination degree of the tube for specific examinations.

"INTELLIGENT"

LCD touch screen panel on the stand column can easily achieve adjusting parameter near the table.



High-performance System Satisfy User Experience

COLUMN TUBE-HOLDER

Column shifting on rails equipped with telescopic arm. With easy movements management thanks to "handle".


PERFORMANCES

Auto-Stitching Function, Autotracking function, Full DICOM.



Lucerna DR 65
Floor Mounted

**Lucerna DR 65
Series**



**Excellent
image
quality**



01 Flexibility

Carefully studied strategies allow GMI to offer a complete and flexible range of products able to satisfy any Client's need, even the most specific ones.



02 Innovation

The activity of R & D department, always alert on new technologies and in close cooperation with renown "Centers of Excellence" is the key to innovative and efficient solutions.



03 Quality

GMI is fully compliant with the required business models and quality standards; our target being maximum customer satisfaction at all times.



04 Tradition

The company's research & development center is in Switzerland and we take full advantage of the country's long history in X-ray engineering competence. In addition, our global operations are kept Swiss working culture is central to our way of working throughout our company.



GMI NET-WORK

The Swiss Accuracy

GMI is a registered trademarks of General Medical Italia Ltd., all rights reserved.

Because of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the GMI sales organization worldwide.

Availability and packaging may vary by country and is subject to change without prior notice.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features which do not always have to be present in individual cases.

Specifications and designs are subject to change without notice GMI reserves the right to modify the design, packaging, specifications and options described herein without prior notice.

Note: Any technical data contained in this document may vary within defined tolerances.



General Medical Italia Ltd

ENGLAND

17 Carlisle Street, W1D 3BU London (UK)
Tel. +44(0) 207 283 7489 - www.gmitalia.co.uk

SWITZERLAND

Piazzetta L. Fontana 4 -6850 - Mendrisio (CH)
Tel +41 (0)76 534 0461 - www.gmint.ch

ITALY

Via Enrico Fermi 26 – 24050 Grassobbio (BG)
Tel +39 035 65 948 11 - www.gmint.ch

