

Lucerna DR Series

Floor & Ceiling mounted DR



*Lucerna DR65
Ceiling*



The new X-RAY generation
Digital Radiography Solution

Lucerna DR65C

State-of-the-art ceiling suspension system

The **Lucerna DR65 Ceiling** combines efficient space use, a digital detector for rapid image capture and reduced radiation, and an advanced software for enhanced image analysis.

This versatile system delivers superior image quality, increased efficiency, and minimized radiation exposure, adapting seamlessly to diverse clinical demands.

It boasts automatic positioning features, a dual FPD configuration, fully automatic stitching, and much more, with unmatched efficiency and accuracy.



Full-Lenght elevating patient bed

Thanks to the motorized FPD, and the ceiling system autotracking features, the **Lucerna DR65C** highly reduces the need for patient repositioning

Tilting Bucky

Both manual and motorized control, allows the user to fulfill their diagnostic needs



HD Tube Touch Display

System features and Patient data are always at a glance



DXRAYDIAGNOST WORKSTATION

- Management of the entire workflow
- Complete control of digital images
- Outstanding image quality
- Automatic APR setting
- Fully automatic Stitching
- DICOM 3.0

**Lucerna DR65
Ceiling**

Lucerna DR65
Floor Mounted




The new X-RAY generation
Digital Radiography Solution

Lucerna DR65F

Multi-purpose Digital X-Ray System

The **Lucerna DR65 Floor Mounted** is the perfect choice for general X-Ray imaging, where efficiency, capabilities, and high resolution low dose exams are a must.

The **Lucerna DR65F** smart design and ease of use guarantees excellent user interaction and patient comfort; integrated touch controls allow for easy system positioning, and the on board touchscreen display makes every information and setting always available.



Full Body Stitching

The Assisted Stitching Function delivers excellent results, by guiding and assisting the user in every step of the procedure

Everything in reach




Thanks to the 10.1" touchscreen display, information about the patient and the system are always in reach and in sight.

- Patient data
- APR and exam details
- Configurable exposure data
- AEC settings
- Dedicated Stitching interface
- System positioning and status
- Autotracking function



***Lucerna DR65
Floor Mounted***

**Lucerna DR65
Series**



**Excellent
image
quality**





01 Flexibility

Carefully studied strategies allow GMI to offer a complete and flexible range of products able to satisfy any Client's need, even the most specific ones.



02 Innovation

The activity of R & D department, always alert on new technologies and in close cooperation with renown "Centers of Excellence" is the key to innovative and efficient solutions.



03 Quality

GMI is fully compliant with the required business models and quality standards; our target being maximum customer satisfaction at all times.



04 Tradition

The company's research & development center is in Switzerland and we take full advantage of the country's long history in X-ray engineering competence. In addition, our global operations are kept Swiss working culture is central to our way of working throughout our company.



The Swiss Accuracy

The Swiss Accuracy

GMI is a registered trademarks of General Medical Italia Ltd., all rights reserved.

Because of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the GMI sales organization worldwide.

Availability and packaging may vary by country and is subject to change without prior notice.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features which do not always have to be present in individual cases.

Specifications and designs are subject to change without notice GMI reserves the right to modify the design, packaging, specifications and options described herein without prior notice.

Note: Any technical data contained in this document may vary within defined tolerances.



General Medical Italia Ltd

ENGLAND

17 Carlisle Street, W1D 3BU London (UK)
Tel. +44(0) 207 283 7489 - www.gmitalia.co.uk

SWITZERLAND

Piazzetta L. Fontana 4 -6850 - Mendrisio (CH)
Tel +41 (0)76 534 0461 - www.gmint.ch

ITALY

Via Enrico Fermi 26 – 24050 Grassobbio (BG)
Tel +39 035 65 948 11 - www.gmint.ch

